

Certificate of Conformity

VWR Catalogue Number	612-9351
Description	Disposable Inoculating loop, 1 µl , Rigid, Dark Green
Lot Number	2002950
Expiry Date (yyyy-mm-dd)	2023-03-31
Country of Origin	Manufactured in Italy
Date of Issue (yyyy-mm-dd)	2023-01-26

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Polystyrene.
Sterilization	Product labelled as sterile is gamma irradiated. Sterilisation procedure is validated in conformity with BS EN ISO 11137-1:2015 and BS EN ISO 11137-2:2015 to guarantee SAL 10 ⁻⁶ .
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.

Non-Pyrogenic Statement	Not applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-06-10 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	No special storage condition required.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature