

Certificate of Conformity

VWR Catalogue Number	115-2831
Description	VWR sterile polyester wiper
Lot Number	A09-0065-211220
Date of Manufacture (yyyy-mm)	2021-12
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2022-02-16

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E and analysis in accordance with IEST-RP-CC025.

Test Item	Spec	Test Result	Conclusion
LPC($\geq 0.5\mu\text{m}$)	$<550 \text{ counts/cm}^2$	471	Pass
APC($\geq 0.5\mu\text{m}$)	$<2050 \text{ counts/m}^2$	41	Pass
NVR(DIW)	$<0.09 \text{ mg/g}$	0.0024	Pass
NVR(DIW)	$<0.02 \text{ g/m}^2$	0.0176	Pass
NVR(IPA)	$<0.2 \text{ mg/g}$	0.1024	Pass
NVR(IPA)	$<0.04 \text{ g/m}^2$	0.0139	Pass
Basic Weight	$135 \pm 8 \text{ g/m}^2$	135.09	Pass
Sorbent Capacity	$\geq 325 \text{ mL/m}^2$	329	Pass
Sorbency Rate(DIW)	$\leq 1 \text{ s}$	0.24	Pass
Silicone	Not Detected	Not Detected	Pass
Amide	Not Detected	Not Detected	Pass
PAEs	Not Detected	Not Detected	Pass

Test Item	Spec	Test Result	Conclusion
Sodium	<1.05 ppm	0.18	Pass
Chloride	<0.9 ppm	0.06	Pass

Product Specifications

Material	100% Polyester.
Sterilization	Gamma-ray irradiated according to ISO 11137:2006, Dose set meet SAL 10 ⁻⁶ with the dose range Min: 29 kGy; Max: 38 kGy.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The above product is latex free.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Room temperature, dry and shady.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.