

Certificate of Conformity

VWR Catalogue Number	114-9324
Description	VWR® Cleanroom Coverall With Hood Sterile, White, S
Lot Number	230310
Date of Manufacture (yyyy-mm-dd)	2023-03-10
Date of Sterilization(yyyy-mm-dd)	2023-03-25
Sterilization Expiry date (yyyy-mm-dd)	2028-03-25
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2023-04-28

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E.

	Test Method	MD	CD
Break Tensile properties(N)	ASTM D5034	137	86
Elongation rate [%]	ASTM D5034	109	112
Tear prosperities	ASTM D5034	106	78

Product Specifications

Material	60g SMS (Spunbonded Polypropylene fabric /Meltblown /Spunbonded Polypropylene fabric)
Sterilization	Products are Gamma irradiated as per ISO 11137:2006 and meets SAL 10 ⁻⁶ with a specified dose of 26.5 – 34.7 kGY.
ATP Assay	Not Applicable.

DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	0- 50 °C, 45% - 65% humidity, Shady
Shelf Life	5 years.