

TECHNICAL DATA SHEET

Test tube PP 5 ml 12,5 x 90 mm

Article no.	3199
Capacity ml	5 ml
Material	Polypropylene (PP)
Temperature range	-196°C to +121°C
Autoclavable	Tubes can be autoclaved in upright position with caps loosened at +121°C, for up to 30 minutes
Colourp	Natural
Diameter	12,5 mm
Height	90 mm
Closure	Screw cap, internal threaded
Tube bottom	round
Sterile	no
Standards and conformity	ISO 2859-1: Sampling and inspection procedures FDA: Resin conforms to FDA 21 CFR 177.1520 USP: Resin conforms to USP Class VI CONEG / RoHS: Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals REACH (SVHC): Plastic is in conformity to REACH standards LATEX: Material is Latex Free BSE / TSE: Material is BSE / TSE Free
Storage conditions	Store at room temperature in normal warehouse conditions - Avoid temperature variations and humidity - Protect from any possible contamination - Protect from any damage to the packaging
Recommended use	- Verify proper cap closure when using biohazard material and / or chemical reagents - Follow chemical resistance chart recommendations - For use in automated equipment, follow the equipment manufacturer's instructions - Should be used only in the gas phase of Liquid Nitrogen