

Certificate of Conformance**MINNCARE™ MINIMUM REQUIRED
CONCENTRATION (MRC)
TEST STRIPS****Part Number:** 78339-000**Lot Number:** 580047**Expiration Date:** 2024-06-12

I certify that the test strips indicate a positive reaction with a solution mixed to a concentration of 475 ± 30 ppm Peracetic acid.

I certify that the test strips did not indicate a positive reaction with a solution mixed to a concentration of 350 ± 30 ppm Peracetic acid.

The statistical analysis and approval test procedure is on file and available for review.

This lot has been reviewed and meets specified requirements.


Quality Analyst
Analysis Date
Quality Assurance
Release Date

Medivators Inc.
14605 28th Avenue North
Minneapolis, MN 55447 U.S.A.
Telephone: 763-553-3300
Fax: 763-553-3387
www.medivators.com