

Material	20839.930
Material description	2-Propanol
Grade	GPR RECTAPUR
CAS Number	67-63-0
Molecular formula	$(\text{CH}_3)_2\text{CHOH}$
Molecular mass	60.10

Characteristics	Specifications
Assay	$\geq 99.0 \%$
Acidity	Passes test
Acidity or alkalinity	Passes test
Benzene and related substances	Passes test
IR Spectrum	Passes test
Peroxides	Passes test
Boiling point	81.5 - 82.6 °C
Density (20/20)	0.785 - 0.789
Density (25/25)	0.783 - 0.790
Refractive index (20°C)	1.376 - 1.378
Evaporation residue	$\leq 20 \text{ ppm}$
Water	$\leq 0.5 \%$

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France

Analysis at raw material release

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.