

CERTIFICATE OF STERILITY AND COMPLIANCE

Method of Sterilization: EO
Sterilization Date: 2025-05, 2025-06
PO No.: 9092
No. of Cartons: 1380 cases
Product Name: DermAssist Latex Surgical Gloves

Item Nbr	Lot Nbr	Qty	EO Concentration	Exposure Time	Temp	Humidity
133550	43-2506350C1	15	750±50 mg/lit	5 Hrs	55±5°C	30%(Min
133600	43-2506350C1	80				
133650	43-2506350C1	360				
133700	43-2506350C1	270				
133750	43-2506350C1	350				
133800	43-2506350C1	205				
133850	43-2506350C1	80				
133900	43-2506350C1	20				

The product described above was manufactured to the following ASTM standards:

ASTM D3577 Standard Specification for Rubber Surgical Gloves

ASTM D5151 Standard Test Method for Detection of Holes in Medical Gloves

ASTM D5712 Standard Test Method for Analysis of Aqueous Extractable Protein
in Natural Rubber and Its Products Using the Modified Lowry
Method

ASTM F1671 Standard Test Method for Resistance of Materials Used in
Protective Clothing to Penetration by Blood-borne Pathogens