
Date of Certificate: 10/22/2025

Catalog Number: 36944-0010
Description: STRLWR, SMPL SPOON, 1 TSP, PS, 10/CASE
Lot Number: J399473
Date of Manufacture: 9/26/2025
Country of Origin: US

Material: High Impact Polystyrene (HIPS)
We confirm that the material manufactured with HIPS conforms to:
21 CFR 177.1640
21 CFR 178.3297

Gamma Irradiation: Irradiation exposure dose (avg.) = 23.5 min - 50.0 max kGy
ANSI/AAMI/ISO 11137 (latest revision), Method 1
No testing is done for Pyrogens / Endotoxins.

CoSP Report Number: 3847101

Product Shelf Life: N/A - Product should remain in its irradiated state as long as outer packaging remains completely intact.

BSE/TSE Statement: Based on the information provided to us, the HIPS resin used in the manufacture of the above product is BSE/TSE free.
During manufacture, this product does not come into contact with any products or processes that may pose a BSE/TSE risk.

Other:

Bel-Art Products hereby certifies that the above listed product has been manufactured and inspected in accordance with documented procedures and specifications, and that it conforms to all applicable requirements.

Dejia Savage
Dejia Savage
Quality Inspector

Certificate of Processing

STERIGENICS 75 Tilbury Road Salem NJ 08079
TEL 856 935-1400 FAX 856 935-0542 www.sterigenics.com

R55480102

10/21/2025 15:40:54 GMT
Page 1 of 1

Customer Name:	Bel-Art Products-NJ	Processing Facility:	Salem	Work Order #	3847101
P.O.#	MD77524			Sales Order #	3911932
23.5-50.0 kGy		BelArt, Gamma Treatment	Irradiation Date/Time:	10/18/2025 14:07:00 GMT	

SO Line #	Qty	UOM	Customer Item Number	Customer Item Description	Customer Lot Number	Customer Load Number
101.000	96	CA	369460000(30)	SAMPLING SPOON, 3Tsp	J397449	NA
	Dose Map		125_10745	TE4-0C3		
102.000	100	CA	369062008(20)	STERILEWARE,SCOOP,RED,PS,8 OZ	J398573	NA
	Dose Map		125_10727-1	0C3 - 0E3x1.085		
103.000	30	CA	369100000(30)	STERILEWARE STERILE SCOOP AND	J396501	NA
	Dose Map		125_10746	0C3 - 0E3 x 1.086		
104.000	120	CA	369042004(20)	STERILEWARE,SCOOP,RED,PS,4 OZ	J397504	NA
	Dose Map		125_10885	0C3 - 0E3 x 1.14		
105.000	160	CA	369460010	Tablespoon 1.0	J396423	NA
	Dose Map		125_10778	0C3 - 0E3x1.047		
106.000	60	CA	369440010(30)	STRLWR,SMPL SPOON,1 TSP,PS,10/	J399473	NA
	Dose Map		125_10877	0C3 - 0E3		
	566	CA	Total			

Quality Test Summary

Op#	Quality Test Description	Minimum Spec	Maximum Spec	Result	Pass/Fail	-----Signed By ----- User	Date /Time
450.00	Maximum Dose	23.5 kGy	50.0 kGy	41.4 KGY	Pass	JWIDGER	10/21/2025 09:47:16 GMT
		Reason Code Test				JAMES WIDGER	
450.00	Minimum Dose	23.5 kGy	50.0 kGy	25.7 KGY	Pass	JWIDGER	10/21/2025 09:47:33 GMT
		Reason Code Test				JAMES WIDGER	

Sterigenics certifies that the materials listed above (as described by the Manufacturer) received the indicated doses within the precision and accuracy of the dosimetry system employed.

Electronically Signed By: LESLIE GLESSNER
Reason: Work Order Completions

Date: 10/21/2025 15:40:17 GMT