

Certificate of Quality

VWR NA Catalogue Number	89221-592
VWR EU Catalogue Number	215-3406S
Description	High Density Polyethylene Narrow Mouth Bottle 250 ml BULK
Lot Number	VDF18J24
Date of Manufacture (yyyy-mm)	2024-10
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2025-04-10

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Test Name	Test Reference	Accuracy	Result
Leak Test	QWI-QC-03	General Insp. Level-1, AQL -1.5	Leak proof
Dimensional	QWI-QC-23	General Insp. Level-1, AQL -1.5	Within tolerance

Product Specifications

Material	Bottle: High Density Polyethylene Closure: Polypropylene Above material is 100% Virgin, Compliance with FDA 21 CFR, USP Class VI passed.
Temperature Resistance	-100 to +120 °C
Sterilization	Not Applicable.
ATP Assay	Not Applicable.

DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Testing is conducted to quality all material resins using ISO 10993 standards for cytotoxicity and have been shown to be non-toxic.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Ambient temperature.
Shelf Life	Not applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.