

Product Name:	VWR-HCLS-12
Description:	VWR Catalog No: 76514-260 VWR Standard Series Glass Door Laboratory Refrigerator 12 Cu. Ft. capacity One swing glass door, self-closing, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range: 1°C – 10°C Keyed door lock Heavy duty door handle Magnetic door gasket for positive seal Three adjustable shelves (adjustable in ½" increments) LED interior lights are safety shielded and switch controlled Swivel Casters (3 ½") Probe access port (3/8") Forced air draft circulation White powder-coated steel exterior and interior High density urethane foam cabinet insulation Hydrocarbon, natural refrigerant (R290) Two year parts and labor warranty, plus an additional three year compressor parts warranty Exterior dimensions: 25" W x 29 ¾" D x 62 ¼" H without casters; 65 ¾" H with casters Cycle defrost 115V, 60 Hz, 3 Amps, 1/5 HP ETL listed Shipping weight: 258 lbs.