

|                      |                              |
|----------------------|------------------------------|
| Material             | 20839.366                    |
| Material description | 2-Propanol                   |
| Grade                | GPR RECTAPUR                 |
| CAS Number           | 67-63-0                      |
| Molecular formula    | $(\text{CH}_3)_2\text{CHOH}$ |
| Molecular mass       | 60.10                        |

| Characteristics                | Specifications        |
|--------------------------------|-----------------------|
| Assay                          | $\geq 99.0 \%$        |
| Acidity                        | Passes test           |
| Acidity or alkalinity          | Passes test           |
| Benzene and related substances | Passes test           |
| IR Spectrum                    | Passes test           |
| Peroxides                      | Passes test           |
| Boiling point                  | 81.5 - 82.6 °C        |
| Density (20/20)                | 0.785 - 0.789         |
| Density (25/25)                | 0.783 - 0.790         |
| Refractive index (20°C)        | 1.376 - 1.378         |
| Evaporation residue            | $\leq 20 \text{ ppm}$ |
| Water                          | $\leq 0.5 \%$         |

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare  
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250  
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.