

Certificate of Conformity

VWR Catalogue Number	See List
Description	Test tube racks
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2021-12-13

Catalogue Number	No. of Holes	For tube Ø [mm]	Colour	W x D x H [mm]
212-0396	9 (3×3)	30	Blue	110.6 × 110.6 × 84.2
212-0397	9 (3×3)	30	Green	110.6 × 110.6 × 84.2
212-0398	9 (3×3)	30	Red	110.6 × 110.6 × 84.2
212-0399	9 (3×3)	30	White	110.6 × 110.6 × 84.2
212-0394	16 (4×4)	25	Blue	123.9 × 123.9 × 74
212-0395	16 (4×4)	25	White	123.9 × 123.9 × 74
212-0390	20 (4×5)	20	Blue	130 × 105 × 84
212-0391	20 (4×5)	20	Green	130 × 105 × 84
212-0392	20 (4×5)	20	Red	130 × 105 × 84
212-0393	20 (4×5)	20	White	130 × 105 × 84
212-0425	24 (3×8)	30	Blue	110 × 285.5 × 86.7
212-0426	24 (3×8)	30	Green	110 × 285.5 × 86.7
212-0427	24 (3×8)	30	Orange	110 × 285.5 × 86.7
212-0428	24 (3×8)	30	Red	110 × 285.5 × 86.7
212-0429	24 (3×8)	30	White	110 × 285.5 × 86.7
212-0430	24 (3×8)	30	Yellow	110 × 285.5 × 86.7
212-0382	36 (6×6)	13	Blue	106 × 106 × 57
212-0383	36 (6×6)	13	Green	106 × 106 × 57
212-0384	36 (6×6)	13	Red	106 × 106 × 57
212-0385	36 (6×6)	13	White	106 × 106 × 57
212-0386	36 (6×6)	16	Blue	131 × 131 × 70

Catalogue Number	No. of Holes	For tube Ø [mm]	Colour	W x D x H [mm]
212-0387	36 (6×6)	16	Green	131 × 131 × 70
212-0388	36 (6×6)	16	Red	131 × 131 × 70
212-0389	36 (6×6)	16	White	131 × 131 × 70
212-0413	40 (4×10)	20	Blue	104 × 253.6 × 86.6
212-0414	40 (4×10)	20	Green	104 × 253.6 × 86.6
212-0415	40 (4×10)	20	Orange	104 × 253.6 × 86.6
212-0416	40 (4×10)	20	Red	104 × 253.6 × 86.6
212-0417	40 (4×10)	20	White	104 × 253.6 × 86.6
212-0418	40 (4×10)	20	Yellow	104 × 253.6 × 86.6
212-0419	40 (4×10)	25	Blue	123 × 301 × 94.5
212-0420	40 (4×10)	25	Green	123 × 301 × 94.5
212-0421	40 (4×10)	25	Orange	123 × 301 × 94.5
212-0422	40 (4×10)	25	Red	123 × 301 × 94.5
212-0423	40 (4×10)	25	White	123 × 301 × 94.5
212-0424	40 (4×10)	25	Yellow	123 × 301 × 94.5
212-0400	72 (6×12)	13	Blue	103.5 × 202.5 × 60
212-0401	72 (6×12)	13	Green	103.5 × 202.5 × 60
212-0402	72 (6×12)	13	Magenta	103.5 × 202.5 × 60
212-0403	72 (6×12)	13	Orange	103.5 × 202.5 × 60
212-0404	72 (6×12)	13	Red	103.5 × 202.5 × 60
212-0405	72 (6×12)	13	White	103.5 × 202.5 × 60
212-0406	72 (6×12)	13	Yellow	103.5 × 202.5 × 60
212-0407	72 (6×12)	16	Blue	127.3 × 250.6 × 73
212-0408	72 (6×12)	16	Green	127.3 × 250.6 × 73
212-0409	72 (6×12)	16	Orange	127.3 × 250.6 × 73
212-0410	72 (6×12)	16	Red	127.3 × 250.6 × 73
212-0411	72 (6×12)	16	White	127.3 × 250.6 × 73
212-0412	72 (6×12)	16	Yellow	127.3 × 250.6 × 73

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** and **ISO 13485:2016** standards. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Test Name	Test Reference	Accuracy	Result
Autoclaving	QWI-QC-07	General Insp. Level-1, AQL -1.5	PASS
Dimensional	QWI-QC-23	General Insp. Level-1, AQL -1.5	Within Tolerance

Product Specifications

Material	30% Glass filled Nylon with color masterbatch.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	Intentional addition of raw materials made from substances of animal origin, but comply with European Regulation (EC) N° 1774/2002 and Directive 2000/6/EC. Also the requirements of paragraph 6.4 Tallow derivatives of EMEA/410/01 rev.3 1 July 2011 are met.
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.

REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Ambient temperature.
Shelf Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.