



Product Specification Sheet

VWR CLEANROOM CONSTRUCTION TAPE CTP™

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Width	Pk
817-0062	25.40 mm	36 Roll
817-0077	50.80 mm	1 Roll

Other Specification:

Material: **Polyethylene with synthetic rubber adhesive**

Backing: **White polyethylene**

Adhesive: **Synthetic rubber**

Adhesion to steel: **66 N/100 mm**

Tensile strength: **315 N/100 mm**

Total thickness: **0.18 mm**

Elongation at break: **80%**

Optimum temperature range: **18 to 77 °C (65 to 170°F)**

Abrasion resistance: **Good**

Length: **32.92 m**

VWR International bvba/sprl, Researchpark Haasrode 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>



Version: 2.0



Product Specification Sheet

VWR CLEANROOM CONSTRUCTION TAPE CTP™

Other Characteristics:

- Clean construction
- Easily bands different mediums
- Serrated edge for easy tear
- Aggressive adhesive

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use
- If there is any doubt ask us.

Date of issue: **October 2012**

VWR International bvba/sprl, Researchpark Haasrode 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>



Version: 2.0