

Certificate of Conformity

VWR Catalogue Number	129-0296
Description	Zip-Seal Bag, Single Pocket, with 150 × 100 mm, Dimension Non-Sterile
Lot Number	030323
Date of Manufacture (yyyy-mm-dd)	2023-03-29
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2023-11-09

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E.

Weld strength : $\geq 15\text{N}/25\text{ mm}$.
Tensile Strength : (Vertical $>10\text{MP}/25\text{ mm}$ / Horizontal $> 20\text{MP}/25\text{ mm}$.
Specific Elongation : (Vertical $\geq 450\%$ / Horizontal $\geq 250\%$.

Product Specifications

Material Low Density Polyethylene (LDPE) C150Y + Linear Low-Density Polyethylene (LDPE) 218W (virgin).

Monomers and additives used to produce the products are listed in List of Authorised Substances of Regulation EU/10/2011.

The material grades conform to EU Directive 1935/2004 on materials and articles intended to come into contact with food.

Sterilization Non-Sterile.

ATP Assay Not applicable.

DNase & RNase Free Not applicable.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Room temperature.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.