

Certificate of Conformity

VWR Catalogue Number	689633-010101
Description	Petri dishes, PS, with 3 vents, 90 mm Diameter, 15 mm Height, Sterile
Lot Number	020524
Date of Manufacture (yyyy-mm-dd)	2024-02-05
Date of Sterilization (yyyy-mm-dd)	2024-02 – 2024-03
Expiry Date (yyyy-mm-dd)	2029-02-05
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2025-08-12

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Polystyrene (PS), Color: Transparent.
Sterilization	Product is ETO treated and they meet a minimum requirement of 10^{-6} Sterility assurance level. Sterilisation procedure is validated in conformity with ISO 11135.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not Applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed).
REACH Statement	Not Applicable.
Storage Conditions	Not Applicable.
Shelf Life	5 years.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Annex

Height of Bottom:	13.20 – 13.80 mm
Diam of top bottom:	87.00 – 88.00 mm
Diam of bottom:	86.50 – 87.10 mm
Height of lid:	7.60 – 8.20 mm
Diam of lid:	90.9 – 91.5 mm