

Certificate of Quality

VWR Catalogue Number	59582-585
Description	VWR® Rubber Stopper Assortment, Black
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2026-04-08

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared the products are manufactured in accordance with their internal quality system. Products are inspected and controlled through whole production processing in accordance with their applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Styrene-Butadiene Rubber, Natural rubber blend (M29).
Sterilization	The product is non-sterile.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	No testing performed.
Non-Pyrogenic Statement	No testing performed.
Latex Statement	The product contains naturally dry latex.
BPA Statement	No testing performed.

DEHP Statement	No testing performed.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP), are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Store in cool, dry climate condition. Store under 75% < rH; 65-75 °F, Keep out of direct UV Light.
Shelf Life	5 years under proper storage conditions.

Disclaimer: VWR'S Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co- manufacturing, or processing activities such as design, packaging, labelling, testing, or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use. This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	April 2026