

Certificate of Conformity

Catalogue Number	49002-31061
Description	CONDUCTIVE TIPS 200 µL FILTER TRAY STERILE CS2304 IVD
Lot no.	1016602/1-1
Country of Origin	Manufactured in Germany
Date of Issue (dd/mm/yyyy)	03/12/2025
Production date (mm/yyyy)	02/2022
Expiry date (mm/yyyy)	02/2027

Quality System Compliance

Products are manufactured in [ISO 13485:2016](#) and [ISO 9001:2015](#) certified plant. Products are inspected and controlled in accordance with the current applicable product specification.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Tips: Polypropylene Filter: Polyethylene C-tray: Polypropylene
Manufacturing Environment	Robotic tips are manufactured under clean room conditions according to DIN EN ISO 14644-1:2016-06.
Sterilization (only for sterile variants)	Products are EB (Electron beam) irradiated and dose released upon ISO 11137-2:2013 recommended practices in effect at the time of validation meet the requirement of 10^{-6} SAL
ATP Assay	Not applicable
PCR inhibitors Free	Product yield equivalent to control < 10 amplifiable targets
Human DNA Free	Tested with detection level < 1 pg
DNase & RNase Free	Tested with detection level < $1 \cdot 10^{-6}$ KunitzU (DNase) and < $1 \cdot 10^{-9}$ KunitzU (RNase)
BSE/TSE	No use of any raw materials produced from or substances derived from animal origin. The manufacturing process of the product does not intentionally use any ingredient of animal origin and no material derived from or exposed to animals affected by.
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	See Endotoxine free

Latex Statement	No intentional use in the formulation of this product.
BPA Statement	No intentional use in the formulation of this product.
DEHP Statement	No intentional use in the formulation of this product.
Endotoxine Free (Only for sterile variants)	Tested with detection level < 0,02 EU/Tip
RoHS	The raw materials used are compliant with the currently valid directive 2011/65/EU and 2002/95/EG (RoHS).
REACH Statement	The product does not contain any "SVHC" ingredients, in accordance with the currently valid regulation (EC) No 1907/2006 on REACH.
Storage Conditions	Keep dry and away from sunlight.
Shelf Life	5 years

Disclaimer: Ritter states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
 This document has been produced electronically and is valid without a signature.



Ritter GmbH
 Kaufbeurer Str. 55
 86830 Schwabmünchen (Germany)
 Phone: +49 8232 5003-45
 E-Mail: medical@ritter-online.de