

Certificate of Conformity

VWR Catalogue Number	See List
Description	Aprons, VWR® Basic, White Colour, 20 µm Thickness
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2021-12-29

Catalogue Number	Dimensions [mm]
113-9094	800 × 1400
113-9093	800 × 1250

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP. Inspection records are viewed and signed off by qualified personnel for product release.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E.

Item		Test Method	Standard
Break Tensile Properties ((N/5cm))	MD	ISO13934-1	>20
	CD	ISO13934-1	>10
Elongation Rate (%)	MD	ISO13934-1	>100
	CD	ISO13934-1	>100

Product Specifications

Material 20 µm Polyethylene film.

Comply with European Commission Directive 1935/2004/EC, Regulation (EU) No.10/2011, Regulation (EU) No.2016/1416 and Regulation (EU) 2017/752

Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Room temperature, dry and shady.
Shelf Life	5 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.