

Certificate of Conformity

VWR Catalogue Number	See List
Description	HDPE Bottle with Plug
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2026-01-06

Catalogue Number	Capacity [mL]	Bottle Shape	Sterility	Cap Colour	Neck Size & finish [mm]	Length [mm]	Width [mm]
215-4497	150	Square	Sterile	Red	45/400	50.5	50.5
215-4498	250	Square	Sterile	Red	45/400	57	57
215-4499	500	Square	Sterile	Red	63/410	76.5	76.5
215-4504	1000	Rectangle	Sterile	Red	63/410	94.9	76.3
215-4505	2000	Rectangle	Sterile	Red	63/410	108	97
215-4511	150	Square	Non-Sterile	Blue	45/400	50.5	50.5
215-4512	250	Square	Non-Sterile	Blue	45/400	57	57
215-4513	500	Square	Non-Sterile	Blue	63/410	76.5	76.5
215-4514	1000	Rectangle	Non-Sterile	Blue	63/410	94.9	76.3
215-4515	2000	Rectangle	Non-Sterile	Blue	63/410	108	97

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015 & ISO 15378:2017**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Test Name	Test Reference	Accuracy	Result
Leak Test	SOP/QA/07A	General Insp. Level-1, AQL -1.5	Leak Proof
Dimensional	SOP/QA/01	General Insp. Level-1, AQL -1.5	Within Tolerance

Product Specifications

Material	Bottle: High Density Polyethylene, 100% Virgin, FDA 21 CFR Compliance, USP Class VI Passed, Colour: Natural. Closure with Plug: Polypropylene, 100% Virgin, FDA 21 CFR Compliance, USP Class VI Passed. Monomers and additives used to produce the products are listed in List of Authorised Substances of Regulation EU/10/2011. The material grades conform to EU Directive 1935/2004 on materials and articles intended to come into contact with food.
Sterilization	Product labelled as sterile is ETO (Ethylene Oxide) irradiated as per ISO 11135:2014 and meets SAL 10 ⁻³ . (Which is not applicable for non-sterile products).
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	This product is free of cytotoxicity effects.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.

RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-06-10 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Nitrosamines	The raw material used to manufacture above product does not contain Nitrosamines.
Melamine	Melamine is not intentionally used and is not expected to be present in the product based on the composition of the raw materials used in the manufacture of the product.
Storage Conditions	Store at 25 + 3 / - 8 °C, Relative Humidity 50 ± 10 % RH.
Shelf Life	2 years.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.1	Created General CoQ using new layout	10-2022
2	2.0	Modified Quality system compliance and disclaimer	07-2025
3	2.1	Shelf life-Removed "Which is not applicable for non-sterile products"	01-2026