

# Certificate of Conformity

|                                         |                                                                                         |
|-----------------------------------------|-----------------------------------------------------------------------------------------|
| <b>VWR Catalogue Number</b>             | 631-0108                                                                                |
| <b>Description</b>                      | Adhesion slides, SuperFrost Plus, VWR® 75 mm length, 25 mm width, 1 mm thickness, white |
| <b>Lot Number</b>                       | 051022-9                                                                                |
| <b>Date of Manufacture (yyyy-mm-dd)</b> | 2022-05-10                                                                              |
| <b>Country of Origin</b>                | Manufactured in USA                                                                     |
| <b>Date of Issue (yyyy-mm-dd)</b>       | 2022-10-10                                                                              |

## Quality System Compliance

Products are manufactured under the **ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

## Product Specifications

|                                |                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Material</b>                | Glass and white SF (Super Frost) Ink.                                                                                                                                                                                                                                                                                                                       |
| <b>Sterilization</b>           | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>ATP Assay</b>               | The product does not contain ATP.                                                                                                                                                                                                                                                                                                                           |
| <b>DNase &amp; RNase Free</b>  | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>BSE/TSE</b>                 | No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE). |
| <b>Cytotoxicity</b>            | The product is not cytotoxic.                                                                                                                                                                                                                                                                                                                               |
| <b>Non-Pyrogenic Statement</b> | Product does not contain Pyrogenic substances.                                                                                                                                                                                                                                                                                                              |

|                           |                                                                                                                                                                                                                                                                                                                                                                      |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Latex Statement</b>    | Not Applicable.                                                                                                                                                                                                                                                                                                                                                      |
| <b>BPA Statement</b>      | Product does not contain BPA.                                                                                                                                                                                                                                                                                                                                        |
| <b>DEHP Statement</b>     | Product does not contain DEHP.                                                                                                                                                                                                                                                                                                                                       |
| <b>RoHS</b>               | No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed. |
| <b>REACH Statement</b>    | No substances as listed on >Substance of Very High Concern< candidate list published 2022-06-10 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.                                                                                                                                                      |
| <b>Storage Conditions</b> | General Laboratory Storage Conditions.                                                                                                                                                                                                                                                                                                                               |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.