

# Certificate of Conformity

|                                   |                                                                                                       |
|-----------------------------------|-------------------------------------------------------------------------------------------------------|
| <b>VWR Catalogue Number</b>       | 111-0615                                                                                              |
| <b>Description</b>                | VWR®, Cleanroom Disposable Goggles, Cleanroom, sterile, anti-fog, anti-scratch, Clear PC, Transparent |
| <b>Country of Origin</b>          | Manufactured in Malaysia                                                                              |
| <b>Date of Issue (yyyy-mm-dd)</b> | 2023-09-15                                                                                            |

## Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 14001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

## Product Specifications

Product complies with regulation (EU) 2016/425 and with the standard EN 166:2001.

|                               |                                                                                                                                                                                                                                                                                                                                                              |
|-------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Material</b>               | Lightweight ultra-soft PVC frame.<br>Toughened polycarbonate lens.<br>Latex-free silicone head band.                                                                                                                                                                                                                                                         |
| <b>Sterilization</b>          | Products are Ethylene oxide sterilized and meets a minimum requirement of 10 <sup>-6</sup> SAL.                                                                                                                                                                                                                                                              |
| <b>ATP Assay</b>              | Not applicable.                                                                                                                                                                                                                                                                                                                                              |
| <b>DNase &amp; RNase Free</b> | Not applicable.                                                                                                                                                                                                                                                                                                                                              |
| <b>BSE/TSE</b>                | No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE). |
| <b>Cytotoxicity</b>           | Not applicable.                                                                                                                                                                                                                                                                                                                                              |

|                                |                                                                                                                                                                                                                 |
|--------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Non-Pyrogenic Statement</b> | Not applicable.                                                                                                                                                                                                 |
| <b>Latex Statement</b>         | This product is not made from natural rubber latex.                                                                                                                                                             |
| <b>BPA Statement</b>           | Bisphenols are not used in the manufacture of the raw material and are not expected to be present.                                                                                                              |
| <b>DEHP Statement</b>          | Not applicable.                                                                                                                                                                                                 |
| <b>RoHS</b>                    | Not applicable.                                                                                                                                                                                                 |
| <b>REACH Statement</b>         | No substances as listed on >Substance of Very High Concern< candidate list published 2023-06-14 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed. |
| <b>Storage Conditions</b>      | Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.                                                                                                                             |
| <b>Shelf Life</b>              | 3 years.                                                                                                                                                                                                        |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.

## Version History

| Sl.no | Version number | Description                          | Month & Year |
|-------|----------------|--------------------------------------|--------------|
| 1     | 1.0            | Created General CoC using new layout | 09-2023      |
|       |                |                                      |              |
|       |                |                                      |              |