

Material	26098CF
Material description	Fluorescein sodium
Grade	TECHNICAL dye water and alcohol soluble
CAS Number	518-47-8

Characteristics	Specifications
Assay	≥ 80.0 %
IR Spectrum	Passes test
UV Spectrum	Passes test
Loss on drying (110°C)(at manufacturing)	≤ 10 %
Absorption maximum (buffer pH 8.0)	490 - 492 nm
Spec.absorptivity(λm;0.005g/l;pH8.0;dry)	≥ 1200

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Jessie Pijck, Head of Laboratory - Haasrode
VWR International bv; Geldenaaksebaan 464; BE-3001
Leuven; Belgium

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.