



AC-18726x	CUSTOMER DRAWING	Page 1 of 1
Revision 02	PP Tube, two-position vent stopper, 14 ml	 greiner bio-one
	Item-No. 187 26x	



NOT TO SCALE



Cat.-No.	Description	A	B	C	D	E	F
187 261	Tube	Ø 15.0	Ø 18.0	Ø 22.2	95.0	95.6	19.0
187 262	Tube	Ø 15.0	Ø 18.0	Ø 22.2	95.0	95.6	19.0

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date 19 January 2004	Date 19 January 2004	Date 20 January 2004	
Date September 13, 2005	Name F. Falter	Name F. Sonntag	Name F. Sonntag	