



Quality Certificate

VWR® DISPOSABLE POLYPROPYLENE SPATULAS

As per the manufacturer, the below product meets the following criteria:

North American Catalog Number:	80081-192
Lot Number:	GD190019
Description:	VWR Spatula PP Disp GRN PK150
Date of Manufacture:	Q2 2020
Country of Origin:	China

FDA Statement: The product is a polypropylene tool for spooning, scooping, scraping and stirring of a variety of samples or chemicals. This tool is for Research use only and is not a medical device. The materials, manufacturing process, standards and procedures used to manufacture the product, meet FDA standards for plastics coming into contact with food as defined by the FDA in 21CFR Part 177 subpart B: Substances for Use as Basic Components of Single and Repeated Use Food Contact Surfaces, part 177.1520

Material used in Manufacturing: 100% Food Grade Polypropylene Homopolymer

BSE/TSE Statement: No raw materials that contain, or are derived from, animals are specified in the manufacture of this product.

Quality Control: Random samples of each batch are reviewed for adherence to size, weight and other manufacturing specifications.

Storage Conditions: If stored in a cool and dry location polypropylene resins should have an unlimited shelf life. Heat and direct sunlight will degrade most plastic polymers. Storage of the above referenced products out of direct sunlight and in normal lab conditions should result in an unlimited shelf life.

Non-Pyrogenic: N/A

DNase & RNase Free: N/A

Cytotoxicity: N/A

Latex Free Statement: N/A

Signed:

Jamie Ethier
VP Global Quality
VWR, Part of Avantor

Date: September 29, 2020