

Certificate of Conformity

Catalogue Number	See list
Description	Robotic tips, non-sterile, blistered or stacked
Country of Origin	Manufactured in Germany
Date of Issue (dd/mm/yyyy)	06/09/2022

Catalogue Number	Item Description	Volume
49000-0006	Conductive tips 300 µL tray	300 µL
49001-0106	Conductive tips 1100 µL filter tray	1100 µL
49002-0006	Conductive tips 200 µL S-tray	200 µL
49002-0106	Conductive tips 200 µL filter S-tray	200 µL
49002-0007	Conductive tips 200 µL C-tray	200 µL
49002-0107	Conductive tips 200 µL filter C-tray	200 µL
49024-0008	Conductive tips 350 µL C-tray stacked	350 µL
49003-0006	Conductive tips 1000 µL S-tray	1000 µL
49003-0106	Conductive tips 1000 µL filter S-tray	1000 µL
49003-0007	Conductive tips 1000 µL C-tray	1000 µL
49003-0107	Conductive tips 1000 µL filter C-tray	1000 µL
49033-0006	Conductive tips 1000 µL wide bore tray	1000 µL
49033-0106	Conductive tips 1000 µL wide bore filter tray	1000 µL
49006-0006	Conductive tips 50 µL S-tray	50 µL
49006-0106	Conductive tips 50 µL filter S-tray	50 µL
49006-0007	Conductive tips 50 µL C-tray	50 µL
49006-0107	Conductive tips 50 µL filter C-tray	50 µL
49008-0003	Conductive tips 300 µL tray	300 µL
49008-0105	Conductive tips 300 µL filter tray	300 µL
49009-0005	Conductive tips 1000 µL tray	1000 µL
49009-0105	Conductive tips 1000 µL filter tray	1000 µL
49010-0003	Conductive tips 50 µL tray	50 µL
49010-0105	Conductive tips 50 µL filter tray	50 µL
49018-0106	Conductive tips 50 µL QG filter tray	50 µL
49209-0105	Conductive tips MP 1000 µL filter tray	1000 µL
49106-0006	Non-conductive tips 50 µL S-tray	50 µL
49102-0006	Non-conductive tips 200 µL S-tray	200 µL
49103-0006	Non-conductive tips 1000 µL S-tray	1000 µL

Quality System Compliance

Products are manufactured in [ISO 13485:2016](#) and [ISO 9001:2015](#) certified plant. Products are inspected and controlled in accordance with the current applicable product specification.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material

Tips: Polypropylene
Filter: Polyethylene

Tray / S-tray - 49000-0006, 49001-0106, 49002-0006, 49002-0106, 49003-0006, 49003-0106, 49033-0006, 49033-0106, 49006-0006, 49006-0106, 49018-0006, 49102-0006, 49103-0006, 49106-0006: Polycarbonate

Tray / C-tray - 49002-0007, 49002-0107, 49024-0008, 49003-0007, 49003-0107, 49006-0007, 49006-0107, 49008-0003, 49008-0105, 49009-0005, 49009-0105, 49010-0003, 49010-0105, 49209-0105: Polypropylene

Manufacturing Environment	Robotic tips are manufactured under clean room conditions according to DIN EN ISO 14644-1:2016-06.
Sterilization	Not applicable
ATP Assay	Tested with detection level $< 1 \times 10^{-12}$
PCR inhibitors Free	Product yield equivalent to control < 10 amplifiable targets
Human DNA Free	Tested with detection level < 1 pg
DNase & RNase Free	Tested with detection level $< 1 \times 10^{-6}$ (DNase) and $< 1 \times 10^{-9}$ (RNase)
Nucleic Acid Free	Not applicable
BSE/TSE	No use of any raw materials produced from, or substances derived from animal origin. The manufacturing process of the product does not intentionally use any ingredient of animal origin and no material derived from or exposed to animals affected by.
Cytotoxicity	Not applicable
Non-Pyrogenic Statement	Not applicable
Latex Statement	Not applicable
BPA Statement	Not applicable
DEHP Statement	Not applicable
RoHS	The raw materials used are compliant with the currently valid directive 2011/65/EU and 2002/95/EG (RoHS).
REACH Statement	The product does not contain any "SVHC" ingredients, in accordance with the currently valid regulation (EC) No 1907/2006 on REACH.
Material Statement	Not applicable
Storage Conditions	Keep dry and away from sunlight
Shelf Life	5 years

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.