



VWR® Chemical Splash Goggles

VWR Chemical Splash Goggles provide optimum comfort and maximum splash protection. These high-quality safety goggles feature a modern, wraparound style; fogless, distortion-free panoramic lenses; indirect venting for liquid splash protection; and internal relief slots to reduce interference with prescription glasses. The wide elastic strap with adjustable buckle attaches to the frame with a swivel clip for a comfortable fit. The nose area is split for greater comfort and the wide facial flange provides a secure seal to most every face. Goggles provide 99.9% UV protection.

Produced under stringent conditions and monitored closely to ensure high levels of quality for every pair. The goggles passed tests according to the American National Standard for Occupational and Educational Personal Eye and Face Protective devices, ANSI Z87.1-2015.

Cat. No. 75982-628	
Features	Indirect Ventilation
	Anti-fog lens coating
	Lightweight PVC frame
	Fully adjustable headband
	Fits over most prescription glasses
	Meets ANSI Z87.1-2015 standards
Material:	
Lens	Polycarbonate (PC), Anti-UV 99.9%
Lens Coating	Anti-Fog
Frame	PVC
Elastic Band	60% Textured Polyester + 40% Rubber Band
Weight and Dimension:	
Total Weight	86.3 g
Lens Thickness	2.0 mm
Lens Base Curves	6
Bridge Width	30 mm
Bridge Height	20 mm
Overall Width	170 mm
Overall Length	75 mm

VWR® Chemical Splash Goggles

Ventilation	Lens Coating	Cat. No.	Unit
Indirect	Anti-Fog, Anti-Scratch	75982-628	Pack of 10, Case of 120