

## Certificate of Quality

**REF**

780210

**LOT**

B2101357

Product description: MASTERBLOCK®, 96 WELL, 1 ML, PP, U-BOTTOM,  
BULK PACKED

Following inspections were carried out and passed:

|              |                                     |        |
|--------------|-------------------------------------|--------|
| raw material | <input checked="" type="checkbox"/> | passed |
| visual       | <input checked="" type="checkbox"/> | passed |
| physical     | <input checked="" type="checkbox"/> | passed |

Free of detectable DNase/RNase/human DNA  
Representative product was tested by Real-Time-PCR.

Detection limits:

DNase  $1.3 \times 10^{-8}$  Kunitz units  
RNase  $5.0 \times 10^{-11}$  Kunitz units  
hDNA 5 pg

Free of detectable endotoxins (pyrogens)  
Representative product was tested according to the United States Pharmacopoeia (USP) Chapter <85>  
Detection limit: 0.03 EU/ml

We are certified according to the Quality System Standard ISO 9001 and ISO 13485.

This certificate is valid for the above mentioned lot and has no validity if no lot number is provided.



Certificate creation date: 2021-01-27

[www.gbo.com](http://www.gbo.com)