

Certificate of Analysis

EU Cat. No. : 733-1363
US Cat. No. : 76620-544
Item Number : 5100850-0500
Item description : dNTP Mix 10mM of each dATP, dCTP, dGTP, dTTP
Lot Number : 24L0206
Date of Issue (dd-mm-yyyy) : 02-12-2024
Use by (dd-mm-yyyy) : 31-01-2026
Store at : -20 °C constant temperature freezer. Avoid multiple freeze/thaws.

Characteristics	Specifications	Measured Values
Functional test	dNTP Mix is tested for amplification of a 613 bp fragment from human genomic DNA using Taq DNA polymerase at standard conditions and must yield a single specific product	Passed
Human DNA	Must not be detected by PCR using primers targeting human genomic DNA	Passed
Volume check	Fill volume $\geq$ 0.5 ml per tube	Passed

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.