



Certificate of Conformity

VWR COLD LIGHT SOURCE, VISILIGHT® LED3

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

VWR Catalog no.: 631-1103

Luminous flux (light guide 1 m × Ø 9 mm)	: 650 lm
Colour temperature	: 5.700 K
Light guide access	: Ø 15 mm
Power consumption	: 44 W
Operating Voltage	: 100 - 240 V AC, 50/60 Hz
LED life span	: 50.000 h (according to manufacturer)
Dimensions	: 259 mm × 174 mm × 81 mm
Material	: anodized aluminum
Colour	: titanium and grey
Weight	: apx. 3.7 kg

Other Characteristics:

Superior finishing quality.

Despite the passive cooling system without fan it easily surpasses the brightness of a 150 W halogen cold light source, which makes this product particularly suitable for clean room applications.

Flicker free illumination due to optimized and low-loss electronics

Robust aluminum housing in modern design

By means of optional adapters suitable for various types of light guides

Vertical and horizontal positioning, stackable

Comes with a desktop power supply

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **July 2015**