



10 X CONCENTRATE DILUENT BUFFER FOR DAB MEMBRANE SUBSTRATE (DAB-10) - DAB-B50

Code: DAB-B50

Size: 50 mL

Product Description: 10 X CONCENTRATE DILUENT BUFFER FOR DAB MEMBRANE SUBSTRATE (DAB-10) - DAB-B50

PhysicalState: Liquid - clear

Label	Unconjugated
Stabilizer	Proprietary 10X TRIS buffer with hydrogen peroxide
Storage Condition	Store container at 4° C prior to opening. Protect from moisture and light. No special shipping conditions or precautions are required.
Application Note	Recommended to be used with DAB-10.
Purity And Specificity	pH: 6.0 +/- 0.1 Stability @ 37°C: PASS Stability @ 18° to 26°C: PASS Stability @ 4°C: PASS QC Raw Material: PASS Absorbance check of final product: PASS Performance Data per Immunohistochemistry: PASS
Immunohistochemistry	1X
WESTERN BLOT	1X
IHC	1X
Expiration	Expiration date is one (1) year from date of opening.

Related Products

610-4302	Anti-MOUSE IgG (H&L) (RABBIT) Antibody Peroxidase Conjugated - 610-4302
611-1302	Anti-RABBIT IgG (H&L) (GOAT) Antibody Peroxidase Conjugated - 611-1302
B304	NORMAL GOAT SERUM (NGS) - B304
BSA-50	BOVINE SERUM ALBUMIN - Fraction V (Immunoglobulin and Protease Free) - BSA-50

Related Links

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 326, Gilbertsville, Pennsylvania, USA.