

PRODUCT DATA SHEET

Discofix®-2 blue (4099010)

Stopcock systems for infusion therapy and monitoring

Attribute Name	Attribute Value
Brand name	Discofix®
Product name (HC)	Discofix®
Product type	Stopcock
Type of stopcock	2-way stopcock
Color	blue
Lipid resistant	Yes
Drug resistant	No
Number of Safeflow Valves	-
max. pressure	2.0 bar
Priming volume	0.200 ml
Rotating adapter	No
Type of packaging	Peel pack
Product (unpack.) contains PVC	No
Product (unpack.) cont. Latex	No

Attribute Name	Attribute Value
Product (unpack.) cont. DEHP	No
Packaging contains PVC	No
Packaging contains Latex	No
Packaging contains DEHP	No
Product contains PVC	No
Product contains Latex	No
Product contains DEHP	No
Product contains Nickel	No
Product contains BPA	Yes
Sterility	sterile
Type of Sterilization	Gamma
Disposable / Reusable	disposable