

Specification

Water bidistilled, sterile

A4042

Physical Description:	Liquid
Product Code:	A4042
Product Name:	Water bidistilled, sterile
Specifications:	Sterility: passes test
WGK:	nwg
Storage:	RT
Molecular Formula:	H ₂ O
M:	18.01 g/mol
CAS:	7732-18-5
EINECS:	231-791-2
CS:	28539010

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 0072250192 • Register court Darmstadt HRB Nr. 7340