

Phenol Liquefied
U.S.P.
Contains 20 ppm Citric Acid as a Stabilizer



Material No.: 0610-06

Revision No.: 0

Product Specification

Meets J.P. Chemical Specifications, Meets U.S.P Requirements, GMP Manufactured Product

Test	Specification
USP - Assay (C_6H_5OH)	$\geq 89.0 \%$
USP - Benzene	≤ 2.0 ppm
USP - Identification A	Passes Test
USP - Identification B	Passes Test
USP - Clarity of Solution and Reaction	Passes Test
USP - Nonvolatile Residue	$\leq 0.05 \%$
Distilling Range - Initial Temperature (For Information Only)	
Distilling Range - Final Temperature	≤ 182.5 °C
Stabilizer	Passes Test
JP - Identification (1)	Passes Test
JP - Identification (2)	Passes Test
JP - Boiling Point	≤ 182 °C
JP - Clarity and Color of Solution and Acidity or Alkalinity	Passes Test
JP - Residue on Evaporation	$\leq 0.05 \%$
JP - Assay	$\geq 88.0 \%$

GMP Manufactured Product

Bulk Pharmaceutical Chemical

CAUTION: For Manufacturing, processing or repackaging

Appearance (colorless to slightly pink liquid)

No Class 2 or 3 solvents are expected to be present. The Class 1 solvent Benzene is used as a raw material in a precursor to phenol. The product is tested for any presence of benzene and the result reported on the Certificate of Analysis.

Elemental Impurities (USP 232) - Information on elemental impurities for this product is available on the associated Product Regulatory Data Sheet and elemental impurity profile report.

Packaging Site: Paris Mfg Ctr & DC


Jamie Ethier
Vice President Global Quality

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700

Avantor Performance Materials, LLC

100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700