



Product Specification Sheet

VWR WAFER BOX SEALING TAPES , 33 m

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Description	Width	Pk
817-0059	Packaging tape	25 mm	36 Roll
817-0063	Packaging tape	50 mm	24 Roll
817-0065	Packaging tape	50 mm	24 Roll

Other Specification:

Material: **Polyethylene**

Backing: **White and pigmented low density polyethylene film**

Adhesive: **Acrylic pressure sensitive**

Adhesion to steel: **62 N/100 mm**

Tensile strength: **192 N/100 mm**

Total thickness: **0.14 mm**

Elongation at break: **400%**

Optimum temperature range:

Operating: -40 to +88 °C (-40 to +190 °F)

Application: 10 to 49 °C (50 to 120 °F)

VWR International bvba/sprl, Researchpark Haasrode 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>



Version: 2.0



Product Specification Sheet

VWR WAFER BOX SEALING TAPES , 33 m

Chemical resistance: **Good**

Abrasion resistance: **Good**

UV resistance: **Good**

Other Characteristics:

- PE film and acrylic adhesive minimise out-gassing and offer good 'clean peel' properties
- Excellent chemical resistance and good conformability
- This low-density PE tape is low in halogens and leachable chlorides.

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use
- If there is any doubt ask us.

Date of issue: **October 2012**

VWR International bvba/sprl, Researchpark Haasrode 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>



Version: 2.0