

Material No.: 76163-598

Material Description: VWR syringe tips, classic style, 12.5mL, 100/pack

Lot No.: 1004566

Date of Manufacturing: 11/30/2021

Expiration Date: 11/30/2026

Compliance and Procedure Result

This manufacturing lot has been sampled and tested in accordance with Standard Operating Procedures and has been released for the following characteristics:

Specifications:

Material: Polypropylene

Manufacturing Environment: Manufactured under clean room conditions according to DIN EN ISO 14644-1:2016-06.

Compliance:

Sterilization: Not applicable

PCR Inhibitor: Product yield equivalent to control <10 amplifiable targets

DNase & RNase Free: Tested with detection level $<1 \times 10^{-6}$ KunitzU (DNase) and $<1 \times 10^{-9}$ KunitzU (RNase)

BSE/TSE: No use of any raw materials produced from or substances derived from animal origin.

The manufacturing process of the product does not intentionally use any ingredient of animal origin and no material derived from or exposed to animals affected by.

Standard: ISO 13485:2016, ISO 9001:2015

Country of Origin: Germany

Signed on behalf of VWR By: Carlin Atkins

Quality Employee

5/16/2023