


CATALOG NO.....		3423-42	
SIZE (NOM VOL. (L)).....		20	
NOMINAL BRIM CAP. (L)....		24.4	
MATERIALS	BOTTLE.....	PC	
	CLOSURE.....	WHITE PP	
	CLOSURE LINER..	SILICONE RUBBER/PP	
NECK FINISH.....		48mm DOUBLE START THR'D	
THREAD (T.P.I.).....		3	
NOM. WEIGHT W/ CLOSURE (gms).....		855	
		INCH	CM
A (NECK I.D. MAX. TOP)...		1.49	3.8
A (NECK I.D. MIN. BTM)...		1.39	3.5
B (HT. W/O CLOSURE).....		19.32	49.1
C (HT. W/ CLOSURE).....		19.42	49.3
D (SQUARE).....		10.05	25.5
W.T. (NOM. WALL THK)....		0.03	0.08

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE MC ²	APP'D PVC 1/9/06	PSS 3423-42	REV. 0106
DESCRIPTION: BIOTAINER, CARBOY, PC, 20L			