

Microbiology & Histology

A word cloud of laboratory equipment and materials. The words are arranged in a horizontal, somewhat rectangular shape. The most prominent words are "Microbiology" and "Histology" in the center. Other words include "Cytology Fixative", "Odor Control", "Clearing Agents", "Needles", "Streakers", "Calibrated", "Non-calibrated", "Formalin", "Neutralizer", "Freezing Media", "Nichrome", "Loops", "Media", "Freezing", and "Cytology Fixative".



460 Glennie Circle
King of Prussia, PA 19406
(800) 332-6647
www.deconlabs.com



Customer Service:
800-932-5000



CitriSolv™ *Solvent and Clearing Agent*

CitriSolv™ is a d-Limonene-based solvent and clearing agent that can be used as a safe alternative to xylene and ethyl acetate. It is soluble with absolute alcohol and will dissolve paraffin, glues, adhesives, and mounting media. CitriSolv provides for better tissue sections – it will not harden or shrink tissue. It dries slowly with no residue and enhances some stains. Use in manual or automatic tissue processor staining procedures in histology/cytology applications that require xylene.

Item	Standard Size	Alternate Size	VWR Catalog #
CitriSolv™	1 G (3.8L)	4 x 1 G (3.8L)	89426-268



CitriSolv™ Hybrid *Solvent and Clearing Agent*

CitriSolv™ Hybrid is a mixture of d-limonene and hydrocarbon material allowing for the reduction of orange aroma and an increase in flash point. Not classified as flammable.

Item	Standard Size	Alternate Size	VWR Catalog #
CitriSolv™ Hybrid	1 G (3.8L)	4 x 1 G (3.8L)	89426-270



HydeAway®

HydeAway® aldehyde decontamination management products effectively convert used aldehydes into non-hazardous material to ensure safety and to facilitate disposal. HydeAway® offers the highest neutralization ratio available and can be used for formaldehydes, formalin and glutaraldehydes. Crystalline HydeAway® turns aldehydes into a solid polymer at a 1:30 ratio and the liquid version turns aldehyde into a

fluid suitable for sewer system disposal at a 1:15 ratio. Tested via EPA methods and proven safe for drain or trash disposal. Treating used aldehyde with Hydeaway™, prior to drain disposal, is better for the environment and makes for a safer workplace.

Item	Standard Size	Alternate Size	VWR Catalog #
Hyde-Away® (Crystals)	30 oz.	12 x 30 oz.	89234-798
Hyde-Away® (Liquid)	1 G (3.8L)	4 x 1 G	89234-796



Flash Freeze™

Flash Freeze™ is moisture-free propelled cold gas used to rapidly cool or freeze samples in clinical or scientific laboratories.

Item	Standard Size	Alternate Size	VWR Catalog #
Flash Freeze™	11 oz.	6 x 11 oz.	89234-810



Inoculating Loops, Streakers and Needles

Made of high quality, twisted nickel chromium wire, LeLoops and LeStreaks are re-usable—making them more economical than disposable plastic loops. Loops slide smoothly across the agar

surface without tearing or gouging. The nichrome wire cools rapidly after flaming, saving time. Loops are available calibrated or non-calibrated. All calibrated loops are subjected to Evan's Blue dye test.

LeLoops come without holders, which are available separately. Holders are made of nickel-plated brass with well insulated handles.

LeStreak and LeStab are permanently attached to lightweight, insulated aluminum holders perfectly balanced for easy manipulation.

LeLoops®				
Item	Loop ID	Std. Size	Alt. Size	VWR Cat #
Calibrated Loop 10 µL		5/pk	5pk/cs	Decon #MP190-25
Calibrated Loop 1 µL		5/pk	5pk/cs	89426-278
Small Loop	2.26mm	25/pk	—	89426-282
Medium Loop	2.91mm	25/pk	—	89426-284
Large Loop	3.57mm	25/pk	—	89426-286
LeStreak®				
Small Streaker Blue	2.26mm	5/pk	—	89426-292
Medium Streaker Green	2.91mm	5/pk	—	89426-294
Large Streaker Red	3.57mm	5/pk	—	89426-296
Calibrated Streaker 1 µL		5/pk	5 x 5/pk	89426-288
Calibrated Streaker 10 µL		5/pk	5 x 5/pk	89426-290
LeStab® Needles				
LeStab Micro Needle		5/pk	—	89410-656
LeLoop®/ Loop Holders				
Brass Loop Holder	6"	1 each	10/pk	89426-274
Brass Loop Holder	8"	1 each	10/pk	89426-276



Odorgone® & The Pill™

Autoclave Deodorizer

Heat and moisture from the autoclave release the contents of the capsule to overcome unpleasant smells generated from heating microbiological waste.

Simply place one capsule in each load to absorb offensive odors leaving the laboratory smelling fresh and clean. Odorgone comes in two pleasant scents; Lemon and Apple. The PILL is fresh scented.

Item	Standard Size	Alternate Size	VWR catalog #
OdorGone—Apple	25 Capsules	8 x 25 Capsules	89426-262
OdorGone—Apple	100 Capsules	6 x 100 Capsules	89426-264
OdorGone—Lemon	25 Capsules	8 x 25 Capsules	89426-258
OdorGone—Lemon	100 Capsules	6 x 100 Capsules	89426-260
The PILL	100 Capsules	20 x 100 Capsules	89426-254
The PILL	500 Capsules	6 x 500 Capsules	Decon #5500



SlideRite™

Cytology Fixative

SlideRite™ is a highly concentrated, water-soluble spray that forms a lasting protective coating over all types of cytological smears. Many cytological specimens are fragile and degenerate quickly, therefore, an immediate fixation of the material to be studied is essential.

Item	Standard Size	Alternate Size	VWR Catalog #
SlideRite	1.5 oz.	12 x 1.5 oz.	89426-266



460 Glennie Circle
King of Prussia, PA 19406
(800) 332-6647
www.deconlabs.com



Customer Service:
800-932-5000