

Material	22081.326
Material description	n-Butyl acetate
Grade	GPR RECTAPUR
CAS Number	123-86-4
Molecular formula	CH ₃ COO(CH ₂) ₃ CH ₃
Molecular mass	116.16

Characteristics	Specifications
Assay	≥ 99 %
IR Spectrum	Passes test
Boiling point	124 - 128 °C
Density (20/4)	0.874 - 0.884

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.