

Certificate of Conformity

VWR Catalogue Number	See List
Description	PURE Image Cleanroom Paper
Country of Origin	Manufactured in South Korea
Date of Issue (yyyy-mm-dd)	2026-07-14

Catalogue Number	Type	Colour	Weight [g/m ²]
115-0149	Non latex paper, A4 size	White	80
115-0150	Non latex paper, A4 size	Blue	80
115-0155	Non latex paper, A4 size	Violet	80
115-0156	Non latex paper, A4 size	White	105
115-1961	Latex-impregnated paper, A2 size	White	104
115-1952	Latex-impregnated paper, A3 size	Blue	85
115-1958	Latex-impregnated paper, A3 size	White	85
115-1954	Latex-impregnated paper, A3 size	Blue	104
115-1959	Latex-impregnated paper, A3 size	White	104
115-0163	Latex-impregnated paper, A4 size	White	80
115-0164	Latex-impregnated paper, A4 size	Blue	80
115-0165	Latex-impregnated paper, A4 size	Yellow	80
115-0166	Latex-impregnated paper, A4 size	Green	80
115-1955	Latex-impregnated paper, A4 size	Pink	85
115-1956	Latex-impregnated paper, A4 size	White	104
115-1957	Latex-impregnated paper, A4 size	White	104
115-1951	Latex-impregnated paper, A4 size	Blue	105
115-1960	Latex-impregnated paper, A5 size	White	104
115-1953	Latex-impregnated paper, A5 size	Blue	105

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015 and ISO 14001:2015**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Specifications		80 g/m²	100 g/m²	Unit	Test method
Components		Pulp (cellulose)		93.10%	
Ingredients		Chemicals		6.90%	
		Water			
		Acryl Resin			
		Antistatic Agent			
		PVA			
Basis Weight		80 ± 3	100 ± 4	g/m2	KSM7013
Caliper		95 ± 5	119 ± 5	µm	KSM7021
Moisture		5.5 ± 0.5	5.5 ± 0.5	%	
pH		7 ± 0.5	7 ± 0.5		
Smoothness	Top side	35 ↑	30 ↑	sec	KSM7028
	Wire side	35 ↑	30 ↑		
Tensile Strength	Machine direction	4.5 ↑	5.0 ↑	kg/15mm	KSM7014
	Cross direction	3.0 ↑	3.5 ↑		
Surface Resistance		<10E11		Ω	Super Megommeter
APC (0.02 CF)	≥ 0.3µm (<3000)	2490		count/ft³	SGS Testing Method Helmke-Yeich Tumble Drum
	≥ 0.5 µm (1000)	501			
	≥ 1.0 µm (-)	103			
Releasable particle and fiber	≥ 0.5 µm	8.1	0.7	Count	IEST-RP-CC003.4
	Nitrate (NO3-)	50 ↓	7.46	ppm(mg/kg)	
	Nitrite (NO2-)	50 ↓	N.D.	ppm(mg/kg)	
	Sulfate (SO42-)	200 ↓	80.4	ppm(mg/kg)	
	Phosphate (PO43-)	50 ↓	0.156	ppm(mg/kg)	
	Floride (F-)	100 ↓	N.D.	ppm(mg/kg)	

Extractable Ions	Chloride (Cl ⁻)	200 ↓	37.3	ppm(mg/kg)	The sample was extracted in DI Water. Analysis was performed by IC
	Bromide (Br ⁻)	50 ↓	1.07	ppm(mg/kg)	
	Potassium (K ⁺)	50 ↓	4.32	ppm(mg/kg)	
	Lithium (Li ⁺)	50 ↓	N.D.	ppm(mg/kg)	
	Sodium (Na ⁺)	300 ↓	65.8	ppm(mg/kg)	
	Calcium (Ca ²⁺)	80 ↓	13.5	ppm(mg/kg)	
	Ammonium	50 ↓	42.0	ppm(mg/kg)	
Heavy Metals	Cadmium (Cd)	0.5 ↓		ppm(mg/kg)	With reference to IEC62321-5 (2013)
	Lead (Pb)	5 ↓		ppm(mg/kg)	
	Mercury (Hg)	2 ↓		ppm(mg/kg)	With reference to IEC62321-4 (2013)
	Hexavalent Chromium Cr (VI)	8 ↓		ppm(mg/kg)	With reference to IEC62321 (2008)
	PBBs	5 ↓		ppm(mg/kg)	With reference to IEC62321-6(2015)
	PBDEs	50 ↓		ppm(mg/kg)	
	BBP(Butyl Benzyl phthalate)	50 ↓		ppm(mg/kg)	With reference to IEC62321-8/CD (2013)
	DBP (Dibutyl phthalate)	50 ↓		ppm(mg/kg)	
	DIBP (Di-isobutyl phthalate)	50 ↓		ppm(mg/kg)	
	DEHP(Di-(2-ethylhexyl phthalate)	50 ↓		ppm(mg/kg)	

Product Specifications

Material	Cleanroom Paper. Suitable for use in Class 10 to Class 1000 (ISO Class 4 to 6) Cleanroom.
Sterilization	Not applicable.
ATP Assay	Not applicable.

DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Store flat in ambient temperature.
Shelf Life	5 Years.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoC using new layout	05-2022
2	2.0	Modified Quality system compliance and disclaimer	09-2025
3	2.1	Specification modified	07-2026