



Product Specification Sheet

VWR PIPETTE TIPS, UNIVERSAL, PP, AUTOCLAVABLE

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Capacity (µl)	Length (mm)	Colour code	Sterile	Packed	Pk
612-5755	2 - 200	50	Yellow	-	Bulk	1.000
612-5756	50 - 1000	70	Blue	-	Bulk	1.000

Material: Polypropylene

Other characteristics:

- Manufactured without lubricants (such as DiHEMDA or oleamide)
- Cadmium-free pigments
- Tips produced in controlled environment and are packaged automatically
- The blue or yellow colour batch does not affect the pipetted media
- Tips are autoclavable at 121°C (2 bar) for 15mins acc. DIN EN 285

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **March 2013**

VWR International bvba/sprl, Researchpark Haasrode 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>

Version: 2.0

