

Material No.: 5621-03
:
Revision No.: 0

Product Specification

Test	Specification
Normality	0.0995 – 0.1005
Sulfate (SO ₄)	≤ 1 ppm
Ammonium (NH ₄)	≤ 3 ppm
Trace Impurities – Heavy Metals (as Pb)	≤ 0.5 ppm
Trace Impurities – Iron (Fe)	≤ 0.1 ppm

For Laboratory, Research, or Manufacturing Use
Standardization at 25°C traceable to NIST Standard Reference Material.

Packaging Site: Paris Mfg Ctr & DC