

# Certificate of Conformity

|                                         |                                                                                |
|-----------------------------------------|--------------------------------------------------------------------------------|
| <b>VWR Catalogue Number</b>             | 201-1399                                                                       |
| <b>Description</b>                      | VWR® Splash Heads, Stutzer, Straight, 29/32 NS socket size, 14/23 NS cone size |
| <b>Lot Number</b>                       | GL0425                                                                         |
| <b>Date of Manufacture (yyyy-mm-dd)</b> | 2025-01-20                                                                     |
| <b>Country of Origin</b>                | Manufactured in India                                                          |
| <b>Date of Issue (yyyy-mm-dd)</b>       | 2025-09-12                                                                     |

## Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015 & ISO 17025:2017**.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

## Product Specifications

The product is in conformance with the requirement of the following standards and regulations:

|                               |                                                                                                                                 |
|-------------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| Technical standards           | NA                                                                                                                              |
| Calibration procedure         | NA                                                                                                                              |
| <b>Material</b>               | The products are made from DIN ISO 3585:1998, TYPE-1, CLASS- A, Boro 3.3 Glass that meet ASTM E-438, TYPE-1, CLASS- A Standard. |
| <b>Sterilization</b>          | Hot air sterilization, in the oven up to 140 °C<br>Steam sterilization, in an autoclave 121 °C/ 20 min<br>134 °C/ 10 min.       |
| <b>ATP Assay</b>              | Not Applicable.                                                                                                                 |
| <b>DNase &amp; RNase Free</b> | Not Applicable.                                                                                                                 |

|                                |                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>BSE/TSE</b>                 | No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE). |
| <b>Cytotoxicity</b>            | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Non-Pyrogenic Statement</b> | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Latex Statement</b>         | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>BPA Statement</b>           | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>DEHP Statement</b>          | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>RoHS</b>                    | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>REACH Statement</b>         | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage Conditions</b>      | Room temperature.                                                                                                                                                                                                                                                                                                                                           |
| <b>Shelf Life</b>              | No predefined shelf life.                                                                                                                                                                                                                                                                                                                                   |

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.