



Technical Datasheet

VWR® BASIC BOUFFANT CAP

CHARACTERISTICS:

- Made of non-woven SPP fabric
- Multiple colours option from white, blue, green, red, yellow to orange.
- Nature Rubber Latex Free
- Particle shedding meet IEST level I requirement
- Comply with European Commission Directive 1935/2004/EC, Regulation (EU) No.10/2011, Regulation (EU) No.2016/1416 and Regulation (EU) 2017/752.


PHYSICAL PROPERTIES

Item		Test method	Standard
Weight(g/m ²)		ISO3801	10 ± 1
Break Tensile Properties (N/5cm)	MD	ISO13934-1	>15
	CD	ISO13934-1	>8
Elongation Rate(%)	MD	ISO13934-1	>30
	CD	ISO13934-1	>30



Technical Datasheet

PRODUCTS AND DIMENSION SPECIFICATION



Unit:CM

ECN		113-1459	113-8762	113-8874	113-8876	113-8875	113-8877	Tolerance
Colour		White	Blue	Green	Red	Yellow	Orange	
Dimension	A	20	20	20	20	20	20	± 1
	A'	53	53	53	53	53	53	± 2
ECN		113-8878	113-8879	113-8880	113-8881	113-8882	113-8883	Tolerance
Colour		White	Blue	Green	Red	Yellow	Orange	
Dimension	A	21	21	21	21	21	21	± 1
	A'	58	58	58	58	58	58	± 2
ECN		113-1460	113-8884	113-8885	113-8886	113-8887	113-8888	Tolerance
Colour		White	Blue	Green	Red	Yellow	Orange	
Dimension	A	22	22	22	22	22	22	± 1
	A'	61	61	61	61	61	61	± 2

Remark: A: Relax, A': Stretch



Technical Datasheet

PACKING:

Package/Quantity	100pcs/bag, 10bags/case
Best Storage Conditions	Temp : 0- 50°C, Humidity: 45-85% without sunshine exposure
Country of Origin	China

QUALITY ASSURANCE CRITERIA:

1. Dimension and Function - Critical (AQL1.5, S-2)

- Product dimension out of spec.
- Fabric characteristics out of spec
- Sonic welding strength less than 15N/5cm

2. Visual Inspection/ Attributes - Major (AQL 2.5, L-II)

- Foreign Material
- Dirt
- Thread ends loose or too long(over 1cm)
- Holes/Tears
- Foul odor
- Discolor
- Wrong collar tag
- Wrong stitch
- Cloth surface applied wrongly
- Less than 3 stitches per centimeter

3. Visual Inspection/Attributes - Minor (AQL 4.0, L-II)

- Sewing sundries
- Elastic band end exposed
- Sewing end falling-off
- Jumping Stitch
- Missing stitch
- Unwanted folds
- Elastic tightness non-consistent

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **November 2015**