



phone 1.847.816.5070
fax 1.847.816.5072
info@heathrowscientific.com
620 Lakeview Parkway
Vernon Hills, IL 60061
www.heathrowscientific.com

TECHNICAL DATA SHEET

Description of Quality Standard: Box is constructed of chipboard cardboard with a white paper covering. Inspected in the assembly facility for incomplete or distorted shape, or other damage or defects, which would compromise proper function. 10 per pack.

Product name:	CARDBOARD CRYOGENIC TUBE STORAGE BOXES AND PARTITIONS
Part number:	HS2860M
Color:	White
Material of construction	Cardboard
Dimensions:	5.7x5.7x4.8"
Packaging:	10-PK

