


Axygen® LTS Compatible Pipet Tips

Effortless ejection, superior sealing,
ergonomic excellence

The new Axxygen LTS Compatible Tips are ergonomically designed to reduce discomfort and optimized to work with Rainin LiteTouch™ pipettors, all with the quality you have come to expect from Corning.

- ▶ **Comfort You Can Feel:** LTS Compatible Tips reduce the force required to seal and eject, enabling a streamlined experience from aspiration to dispensing. Avoid slamming with excessive force on the pipettor for a proper seal, eliminate the struggle to eject tips, and minimize hand and wrist strain after extensive pipetting.
- ▶ **Convenience of Choice:** LTS Compatible Tips are available in your favorite packaging formats, including racked, racked/filtered, tip station, reload stacks, and our new more sustainable solution: Axxygen HybridRack™ Pipet Tips. If you're not using an LTS pipettor, no problem. We provide a complete range of universal pipet tips to choose from.
- ▶ **Quality:** Our United States manufacturing facility is audited by third parties to ensure adherence to regulated manufacturing standards, including ISO 13485 and ISO 14001. Sterile products feature a Sterility Assurance Level (SAL) of 10^{-3} , and all products are made using USP Class VI resin.
- ▶ **Peace of Mind:** Our LTS Compatible pipet tips are RNase-/DNase-free, Human gDNA-free, PCR inhibition-free, and nonpyrogenic so you can pipet stress-free, knowing your samples are safe. Additionally, the tips are autoclavable at 121°C for 15 minutes.



Ordering Information

Products may not be available in all markets.

Axygen® LTS Compatible Pipet Tips

Racked, Filtered, Sterile

VWR Cat. No.	Corning Cat. No.	Description	Qty/Pk	Qty/Cs
77395-566	LTF-20-R-S	20 µL LTS compatible tip, filtered, racked, sterile	960	4,800
77395-568	LTF-200-R-S	200 µL LTS compatible tip, filtered, racked, sterile	960	4,800
77395-570	LTF-1000-R-S	1000 µL LTS compatible tip, filtered, racked, sterile	768	3,840

Racked, Sterile

77395-572	LT-20-R-S	20 µL LTS compatible tip, racked, sterile	960	4,800
77395-574	LT-200-R-S	200 µL LTS compatible tip, racked, sterile	960	4,800
77395-576	LT-1000-R-S	1000 µL LTS compatible tip, racked, sterile	768	3,840

Racked

77395-578	LT-20-R	20 µL LTS compatible tip, racked	960	4,800
77395-580	LT-200-R	200 µL LTS compatible tip, racked	960	4,800
77395-582	LT-1000-R	1000 µL LTS compatible tip, racked	768	3,840

Axygen HybridRack™, Filtered, Sterile

77395-598	LTFER-20-R-S	20 µL LTS compatible tip, filtered, Axxygen HybridRack, sterile	960	4,800
77395-600	LTFER-200-R-S	200 µL LTS compatible tip, filtered, Axxygen HybridRack, sterile	960	4,800
77395-602	LTFER-1000-R-S	1000 µL LTS compatible tip, filtered, Axxygen HybridRack, sterile	960	4,800

Tip Station, Sterile

77395-592	LT-20-TS-S	20 µL LTS compatible tip, tip station, sterile	960	4,800
77395-596	LT-200-TS-S	200 µL LTS compatible tip, tip station, sterile	960	4,800

Tip Station

77395-590	LT-20-TS	20 µL LTS compatible tip, tip station	960	4,800
77395-594	LT-200-TS	200 µL LTS compatible tip, tip station	960	4,800

Reload Stack

77395-584	LT-20-RL	20 µL LTS compatible tip, reload stack	960	3,840
77395-586	LT-200-RL	200 µL LTS compatible tip, reload stack	960	3,840
77395-588	LT-1000-RL	1000 µL LTS compatible tip, reload stack	768	3,072



Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

For additional product or technical information, visit vwr.com/corning, call 1.800.932.5000, or contact your Avantor representative.

avantior™
delivered by **VWR™**

VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2024 Avantor, Inc. All rights reserved.