

|                      |                                                                 |
|----------------------|-----------------------------------------------------------------|
| Material             | 23367.550                                                       |
| Material description | Dichloromethane                                                 |
| Grade                | GPR RECTAPUR Stabilized with about 0,002 % of 2-methyl-2-butene |
| CAS Number           | 75-09-2                                                         |
| Molecular formula    | CH <sub>2</sub> Cl <sub>2</sub>                                 |
| Molecular mass       | 84.93                                                           |

| Characteristics         | Specifications |
|-------------------------|----------------|
| Assay                   | ≥ 99 %         |
| Stabilizer (Isoamylene) | Confirmed      |
| IR Spectrum             | Passes test    |
| Boiling point           | 39.0 - 40.5 °C |
| Acidity                 | ≤ 0.0004 meq/g |
| Density (20/4)          | 1.322 - 1.327  |
| Evaporation residue     | ≤ 10 ppm       |
| Free chlorine           | ≤ 0.3 ppm      |
| Heavy metals (as Pb)    | ≤ 2 ppm        |
| Water                   | ≤ 0.02 %       |
| Cl (Chloride)           | ≤ 3 ppm        |

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare  
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250  
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.