



- Easy and very convenient to use
- Competitive price
- High Quality, Accurate testing result
- ATC Compensation Range(Built-in Automatic Temperature Compensation System): from 10°C to 30°C

Features

1. Automatically measure temperature every 3 seconds, default in centi degree.
2. ATC: 10°C to 30°C (50°F to 86°F).
3. Automatically shutdown after 90 seconds of inactivity.
4. Smaller and lighter, economical version refractometer

Product Code	Product Description	Range	Min. Div.	Accuracy	Remark
RHB-10 ATC (VWRI635-0629)	HANDHELD REFRACTOMETER BRIX 0-10% ATC	Brix:0-10%	Brix: 0.1%	Brix: $\pm 0.10\%$	ATC
RHB-18 ATC (VWRI635-0630)	HANDHELD REFRACTOMETER BRIX 0-18% ATC	Brix:0-18%	Brix: 0.1%	Brix: $\pm 0.10\%$	ATC
RHB-32 ATC (VWRI635-0631)	HANDHELD REFRACTOMETER BRIX 0-32% ATC	Brix:0-32%	Brix: 0.2%	Brix: $\pm 0.20\%$	ATC
RHB-62 ATC (VWRI635-0632)	HANDHELD REFRACTOMETER BRIX 28-62% ATC	Brix:28-62%	Brix: 0.2%	Brix: $\pm 0.20\%$	ATC
RHB-82 ATC (VWRI635-0633)	HANDHELD REFRACTOMETER BRIX 45-82%ATC	Brix:45-82%	Brix: 0.5%	Brix: $\pm 0.50\%$	ATC
RHB-55 ATC (VWRI635-0634)	HANDHELD REFRACTOMETER BRIX:0-55% ATC	Brix:0-55%	Brix: 0.2%	Brix: $\pm 0.20\%$	ATC
RHB-92 ATC (VWRI635-0635)	HANDHELD REFRACTOMETER BRIX 58-92% ATC	Brix:58-92%	Brix: 0.2%	Brix: $\pm 0.20\%$	ATC
RHS-10 ATC (VWRI635-0636)	HANDHELD REFRACTOMETER SALINIT 0-100‰ ATC	Salinity:0-10% Salinity:0-100‰ Specific gravity: 1.000-1.070	Salinity: 0.1% Salinity:1‰ Specific gravity: 0.001	$\pm 0.10\%$ Salinity: ± 0.001 Specific gravity: ± 0.001	ATC
RHS-28 ATC (VWRI635-0637)	HANDHELD REFRACTOMETER SALINITY 0-28%ATC	Salinity:0-28%	Salinity: 0.2%	$\pm 0.20\%$	ATC

RHA-503 ATC (VWRI635-0638)	HANDHELD REFRACTOMETER BATT -50-0° C ATC	EG:-50°C-0°C PG: -50°C-0°C CF: -40°C-0°C BA:1.100-1.400sg	EG:5°C PG: 5°C CF: 5°C BA:0.01sg	EG:±5°C PG: ±5°C CF: ±5°C BA:±0.01sg	ATC
RHB-32S ATC (VWRI635-0639)	HANDHELD REFRACTOMETER ° OE0-140° KMW ATC	0-32% Mas Sacch 0-140°Oe 0-27° KMW Babo	0.2% Mas Sacch 1.0°Oe 0.2° KMW Babo	±0.2% Mas Sacch ±1.0°Oe ±0.2° KMW Babo	ATC



1.800.932.5000 | vwr.com



Prices and product details are current when published; subject to change without notice. Certain products may be limited by federal, state, provincial, or local regulations. VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2017 VWR International, LLC. All rights reserved.