

VWR[®] Rocking Platform Shakers

- 01. FIVE-YEAR WARRANTY
- 02. QUIET OPERATION
- 03. IDEAL FOR WESTERN BLOTTING APPLICATIONS



Produce the precise level of agitation that blotted membranes and gels require

VWR® ROCKING PLATFORM SHAKER

- Easily adjustable tilt range
- 10-120 minute continuous timer
- Safe for coldrooms

An effective speed range and an adjustable rocking motion of $\pm 20^\circ$ provide a gentle “wave” effect that lifts and washes over gels and blotted membranes. Shakers maintain constant speed regardless of load or voltage fluctuation. Choose from rockers with single or double-tiered platforms. Height between platforms on double-tiered models is 9cm (3.5”). Platforms are constructed of stainless steel for durability and easy cleanup. Powerful motor and heavy-duty bearings ensure a long service life. Safe for coldroom use.



Specifications

Ambient Operating Range	-10 to 70°C, up to 85% RH
Maximum Load	5 kg (11 lbs)
Platform Dimensions, W x D	36 x 28 cm (14.2 x 11")
Tilt Range	$\pm 2.5^\circ$ - 20° from horizontal in 2.5° steps
Timer	10 - 120 min, continuous
Speed Range	10 - 120 rpm
Overall Dimensions, Single Platform	36 x 36 x 17 cm (14.2 x 14.2 x 6.7") 5.7 kg (14.3 lbs)
Overall Dimensions, Double Platform	36 x 36 x 26 cm (14.2 x 14.2 x 10.2") 7.8 kg (17.2 lbs)

Description	Cat. No.
VWR Rocking Platform Shaker, Single Platform, 115V	10860-660
VWR Rocking Platform Shaker, Double Platform, 115V	10860-780
Rubber Mat for Platform	75780-066
Additional 2nd Platform	75780-068
Post Kit for 2nd Platform	75780-085