



ZPACK13247037 PLASTIC BUCKET 30 L UN Z



SPECIFICATIONS

Material : HDPE

Exact volume : 37 Liters

Length : 387,00 mm

Width : 286,00 mm

Height : 445,00 mm

Remark : Packaging materials and specifications shown may be modified without preceding notice
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343