

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Pean Hemostatic Forceps
Country of Origin	Manufactured in Pakistan
Date of Issue (yyyy-mm-dd)	2022-01-14

VWR Catalogue Number	Length	Jaw type	Curvature
10806-206	457 mm (18")	Full Serrated Jaw	Straight
10806-210	508 mm (20")	Full Serrated Jaw	Straight
10806-188	184 mm (7.25")	Full Serrated Jaw	Curved
10806-186	184 mm (7.25")	Full Serrated Jaw	Straight
10806-208	457 mm (18")	Full Serrated Jaw	Curved
10806-192	203 mm (8")	Full Serrated Jaw	Curved
10806-190	203 mm (8")	Full Serrated Jaw	Straight
10806-130	610 mm (24")	Full Serrated Jaw	Straight
10806-196	254 mm (10")	Full Serrated Jaw	Curved
10806-128	508 mm (20")	Full Serrated Jaw	Curved
10806-194	254 mm (10")	Full Serrated Jaw	Straight
10806-200	305 mm (12")	Full Serrated Jaw	Curved
10806-198	305 mm (12")	Full Serrated Jaw	Straight
10806-134	610 mm (24")	Full Serrated Jaw	Curved
10806-182	159 mm (6.25")	Smooth Jaws	Straight
10806-204	406 mm (16")	Full Serrated Jaw	Curved
10806-202	406 mm (16")	Full Serrated Jaw	Straight
10806-184	184 mm (7.25")	Smooth Jaws	Straight

Quality System Compliance

Products are manufactured under the **ISO 9001:2008 & ISO13485:2003 & F-899** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	AISI 410 type, Surgical Grade Stainless Steel.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The product is Latex free.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Store in a dry place.
Shelf Life	Not applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.